




SUNDIAL CO. LTD

ΚΑΤΑΛΟΓΟΣ ΠΡΟΙΟΝΤΩΝ
ΕΤΟΙΜΑ ΠΡΟΣ ΧΡΗΣΗ ΘΡΕΠΤΙΚΑ ΥΛΙΚΑ ΜΙΚΡΟΒΙΟΛΟΓΙΑΣ



1. ΕΤΟΙΜΑ ΠΡΟΣ ΧΡΗΣΗ ΘΡΕΠΤΙΚΑ ΥΛΙΚΑ ΜΙΚΡΟΒΙΟΛΟΓΙΑΣ ΓΙΑ ΚΛΙΝΙΚΑ ΕΡΓΑΣΤΗΡΙΑ

Πιστοποίηση & Ποιοτικά Πρότυπα:

Συμμόρφωση με το πρότυπο ISO13485:2016 για Ιατροτεχνολογικά Προϊόντα, απαιτήσεις CE και Ευρωπαϊκού Κανονισμού ΕΕ 2017/746 για IVDR.

DESCRIPTION	REFERENCE	PACKING	MAX. SHELF LIFE
BILE ESCULIN AGAR Selective medium for isolation and identification of Enterococcus	SUN90540	Plates: 10 × 90mm	4 months
BRAIN HEART INFUSION BROTH Medium for cultivation of fastidious aerobic and anaerobic bacteria	SUN15030	Tubes: 20 × 7mL	10 months
BRILLIANT GREEN AGAR Selective medium for <i>Salmonella</i> spp. isolation other than <i>S. typhi</i> .	SUN90550	Plates: 10 × 90mm	3 months
CANDIDA PLUS CHROM AGAR Chromogenic selective medium for the isolation and identification of Candida spp. Including Candida Auris	SUN90210	Plates: 10 × 90mm	3 months
CHOCOLATE MODIFIED w/ ISOVITALEX AGAR Enriched medium for isolating and cultivation of pathogenic microorganisms, in particular for the Neisseria and the Haemophilus species from clinical specimens	SUN90041	Plates: 10 × 90mm	2 months
CLED w/ ANDRADE'S INDICATOR AGAR Medium for total bacterial count of urine	SUN90050	Plates: 10 × 90mm	3 months
CLED WITH ANDRADE'S INDICATOR / MACCONKEY (ENTEROCOCCUS GROWTH) AGAR Differential media for identification, differentiation and confirmation of enteric bacteria from urine specimens. It is appropriate for urine culture like URICULT ^R (dip-slide type).	SUN45420	2-comp. Plates: 10 × 90mm Only on demand Min. quant. 10packs	3 months
COLUMBIA w/ 5% SHEEP BLOOD AGAR General purpose culture medium for the growth and recovery of fastidious microbial species	SUN90470	Plates: 10 × 90mm	2 months
COLUMBIA w/ 5% HORSE BLOOD AGAR General purpose culture medium for the growth and recovery of fastidious microbial species	SUN90020	Plates: 10 × 90mm	2 months
COLUMBIA (w/5% HORSE BLOOD) / CHOCOLATE MODIFIED w/ ISOVITALEX AGAR Selective medium for clinical samples such as throat swabs and generally samples from respiratory system	SUN45231	2-comp. Plates: 10 × 90mm Only on demand Min. quant. 10packs	2 months
COLUMBIA (w/5% HORSE BLOOD) / MAC CONKEY w/ Crystal Violet AGAR Selective media for clinical samples such as vaginal swabs, semen and wounds	SUN45970	2-comp. Plates: 10 × 90mm	2 months
COLUMBIA (w/5% HORSE BLOOD) / MAC CONKEY enterococcus growth AGAR Selective media for clinical samples such as vaginal swabs, semen and wounds	SUN45370	2-comp. Plates: 10 × 90mm	2 months
COLUMBIA (w/5% HORSE BLOOD) / UTI CHROM AGAR Selective media for clinical samples such as vaginal swabs, semen and wounds	SUN45250	2-comp. Plates: 10 × 90mm Only on demand Min. quant. 10packs	2 months
COLUMBIA CNA w/ 5% SHEEP BLOOD AGAR Selective medium for the isolation of Gram(+) bacteria in clinical specimens	SUN90471	Plates: 10 × 90mm Only on demand Min. quant. 5packs	2 months

DERMATOPHYTE TEST AGAR Selective medium for the isolation of pathogenic fungi	SUN15140	Slanted tubes: 10 × 10mL	10 months
	SUN90140	Plates: 10 × 90mm Only on demand Min. quant. 10packs	4 months
DESOXYCHOLATE CITRATE AGAR Selective medium for pathogenic enterobacteria isolation <i>Salmonella</i> and <i>Shigella</i>	SUN90060	Plates: 10 × 90mm	4 months
FASTIDIOUS ANAEROBE AGAR Medium for the cultivation of anaerobic bacteria from clinical and non- clinical specimens	SUN90720	Plates: 10 × 90mm	2 months
GLUCOSE OF MEDIUM Medium recommended for the identification of Gram (-) bacteria based on their ability oxidize or ferment glucose	SUN15611	Tubes: 10 x 5 mL	10 months
	SUN15610	Tubes: 10 x 10 mL Only on demand Min. quant. 10packs	
HEKTOEN ENTERIC AGAR Selective medium for pathogenic enterobacteria isolation <i>Salmonella</i> and <i>Shigella</i>	SUN90290	Plates: 10 × 90mm	3 months
HEKTOEN ENTERIC /S. S. / X.L.D AGAR Triple control for pathogenic enterobacteria isolation <i>Salmonella</i> and <i>Shigella</i>	SUN30270	3-comp. Plates: 10 × 90mm	3 months
MACCONKEY AGAR Selective medium for the isolation of <i>Enterobacteriaceae</i> . Gram positive organisms will not grow on this medium	SUN90100	Plates: 10 × 90mm	4 months
MACCONKEY (ENTEROCOCCUS GROWTH) AGAR Selective medium for the isolation of <i>Enterobacteriaceae</i> including <i>Enterococcus</i> , in biological specimens	SUN90090	Plates: 10 × 90mm	3 months
MACCONKEY w/ crystal violet / S.S. AGAR Selective media for identification, differentiation and confirmation of enteric bacteria from clinical specimens especially fecal specimens	SUN45960	2-comp. Plates: 10 × 90mm	3 months
MACCONKEY (ENTEROCOCCUS GROWTH) / UTI CHROM AGAR Differential media for identification, differentiation and confirmation of enteric bacteria from urine specimens which may contain large number of <i>Proteus</i> species as well as potentially pathogenic gram-positive organisms	SUN45400	2-comp. Plates: 10 × 90mm	4 months
MACCONKEY BROTH Nourishing broth suitable for growing Gram (-) negative and positive for the fermentation of lactose bacteria	SUN15300	Tubes: 20 x 5 mL Only on demand Min. quant. 10packs	10 months
MANNITOL SALT AGAR Selective medium for pathogenic staphylococci isolation	SUN60110	Plates: 10 × 55mm	4 months
	SUN90110	Plates: 10 x 90mm	
MUELLER HINTON II AGAR Recommended for susceptibility test, low concentration of thymine, thymidine and controlled levels of Ca ⁺⁺ and Mg ⁺⁺	SUN90130	Plates: 10 × 90mm	70 days
	SUN12130	Plates: 120x120mm	
MUELLER HINTON w/ sheep blood AGAR Recommended for susceptibility test, for β-hemolytic streptococcus	SUN90132	Plates: 10 × 90mm Only on demand Min. quant. 5packs	2 months
MUELLER HINTON w/ horse blood + 20mg/L β-NAD AGAR Recommended for susceptibility test, for streptococcus spp., haemophilus spp. And other fastidious microbes to antibiotics	SUN90133	Plates: 5 × 90mm Only on demand Min. quant. 5packs	2 months
RAMBACH CHROM AGAR Selective and differential chromogenic culture medium primarily used for the detection and isolation of <i>Salmonella</i> species from complex, multi-microbial samples.	SUN60180	Plates: 10 × 55mm	5 months
RAMBACH CHROM / E. COLI O157 CHROM AGAR Media recommended for the selective isolation of <i>E.coli</i> and enteropathogenic <i>E.coli</i> O157	SUN45410	2-comp. Plates: 10 × 90mm	5 months

SABOURAUD w/ 0,05% CAF AGAR Selective medium for yeasts and moulds isolation	SUN60200 SUN90200 SUN15200	Plates: 10 × 55mm Plates: 10 × 90mm Tubes: 20 × 8mL	4 months 4 months 10 months
SABOURAUD w/ 0,5% CAF AGAR Selective medium for isolation of yeasts and moulds (more selective than Sabouraud Dextrose Agar w/0.05% Chloramphenicol)	SUN60205 SUN90205 SUN15205	Plates: 10 × 55mm Plates: 10 × 90mm Tubes: 20 × 8mL	4 months 4 months 10 months
SABOURAUD CAF + ACTIDIONE AGAR Selective medium for the isolation of fungi, including dermatophytes.	SUN15204	Tubes: 20 x 8 mL	10 months
SALMONELLA-SHIGELLA AGAR Selective medium for <i>Salmonella</i> and <i>Shigella</i> isolation	SUN90190	Plates: 10 × 90mm	3 months
SELENITE BROTH Liquid medium for selective enrichment of <i>Salmonella</i> spp.	SUN70210 SUN15210	Flask: 100mL Tubes: 20 × 7mL	10 months
SCHAEDLER w/ VITAMIN K & 5% SHEEP BLOOD AGAR Medium recommended for the isolation of anaerobic bacteria from clinical specimens	SUN90853	Plates: 10 × 90mm	2 months
SCHAEDLER BROTH w/ VITAMIN K Nourishing broth suitable for growing of demanding aerobic and anaerobic microorganisms	SUN15854	Tubes: 20 x 7 mL	8 months
THIOGLYCOLLATE BROTH w/ vitamin K Medium acc. to Brewer is used for testing sterility of biological products and for isolation of aerobic and anaerobic organisms	SUN15851	Tubes: 10 × 10mL Only on demand Min. quant. 5packs	10 months
TRIPLE SUGAR IRON AGAR Differential medium for the identification of enterobacteria	SUN15700	Tubes: 10 × 10mL	8 months
UTI CHROM AGAR Chromogenic differential medium for identification, differentiation and confirmation of enteric bacteria from specimens such as urine which may contain large number of <i>Proteus</i> species as well as potentially pathogenic gram-positive organisms	SUN90160	Plates: 10 × 90mm	4 months
UREA AGAR Medium for detection of urease production in proteus spp.	SUN15910	Tubes: 10 × 7mL Only on demand Min. quant. 10packs	10 months
YERSINIA C.I.N. AGAR Medium enriched with Cefsulodin, Irgasan and Novobiocin antibiotics for the isolation of <i>Yersinia Enterocolitica</i> in various clinical and non- clinical specimens	SUN90870	Plates: 10 × 90mm	3 months
X. L.D (XYLOSE LYSINE DEOXYCHOLATE) AGAR Differential and selective medium for pathogenic enterobacteria isolation	SUN90280	Plates: 10 × 90mm	3 months
NUTRIENT BROTH suitable for the determination of non fastidious microorganisms	SUN15360	Tubes: 20 x 5 mL	10 months
NUTRIENT AGAR Nutrient rich substrate, suitable for the identification of non fastidious microorganisms	SUN90450	Plates: 10 x 90mm	5 months
THAYER MARTIN AGAR Enriched medium for isolating and cultivation of Neisseria from clinical specimens	SUN90380	Plates: 10 × 90mm	3months

2. ΕΤΟΙΜΑ ΠΡΟΣ ΧΡΗΣΗ ΘΡΕΠΤΙΚΑ ΥΛΙΚΑ ΜΙΚΡΟΒΙΟΛΟΓΙΑΣ ΓΙΑ ΠΕΡΙΒΑΛΛΟΝΤΙΚΑ ΕΡΓΑΣΤΗΡΙΑ, ΝΕΡΟΥ & ΤΡΟΦΙΜΩΝ

Πιστοποίηση & Ποιοτικά Πρότυπα:

Συμμόρφωση με το διεθνές πρότυπο διαχείρισης ποιότητας ISO9001:2015.

DESCRIPTION	REFERENCE	PACKING	MAX. SHELF LIFE
BACILLUS CEREUS AGAR (PREP) Medium for selective isolation and enumeration of <i>Bacillus cereus</i>	SUN90580	10 × 90mm	3 months
BAIRD PARKER AGAR Selective medium for isolation of coagulase positive staphylococci from water and food	SUN60010 SUN90010	10 × 60mm 10 × 90mm	4 months 4 months
BILE ESCULIN AZIDE AGAR Selective medium for group D streptococci isolation	SUN90540 SUN60540	10 × 90mm 10 × 55mm	4 months 4 months
BRAIN HEART INFUSION BROTH Medium for cultivation of fastidious aerobic and anaerobic bacteria	SUN15030	Tubes: 20× 7mL	10 months
BRILLIANT GREEN AGAR (Modified) Selective medium for <i>Salmonella spp.</i> isolation other than <i>S. typhi</i> .	SUN90550	10 × 90mm	4 months
BRILLIANT GREEN BILE BROTH 2% Selective liquid medium for detection and confirmation of coliforms	SUN15600	Tubes: 10× 10mL	8 months
BUFFERED PEPTONE WATER Medium for <i>Salmonella</i> pre-enrichment in food	SUN75460	Bottle: 225mL	8 months
E.C.C CHROM AGAR Chromogenic medium for detection and enumeration of <i>E. coli</i> and other Coliforms in food or water samples	SUN60070	10 × 55mm	4 months
LEGIONELLA BCYE AGAR (w/L-CYSTEINE) Selective medium for isolation of <i>Legionella spp.</i>	SUN60590 SUN90590	10 × 55mm 10 × 90mm	5 months 5 months
LEGIONELLA BCYE AGAR (w/o L-CYSTEINE) Medium for presumptive identification of <i>Legionella spp.</i>	SUN60591 SUN90591	10 × 55mm 10 × 90mm	5 months 5 months
LEGIONELLA GVPC AGAR Medium used in combination with LEGIONELLA BCYE AGAR (w/L-CYSTEINE) and LEGIONELLA BCYE AGAR (w/o L-CYSTEINE) for maintenance of <i>Legionella spp.</i>	SUN60592 SUN90592	10 × 55mm 10 × 90mm	5 months 5 months
LISTERIA CHROM AGAR (ALOA) Chromogenic medium for presumptive identification of <i>Listeria monocytogenes</i> in food	SUN90530	10 × 90mm	5 months
MALT EXTRACT AGAR Medium for the detection, isolation and enumeration of yeasts and moulds	SUN90690	10 × 90mm	3 months
MAXIMUM RECOVERY DILUENT Isotonic diluent for microorganisms' optimal recovery	SUN15390 SUN70390 SUN75390	Tubes: 10 × 9mL Bottle: 100mL Bottle: 225mL	10 months 10 months 10 months
NUTRIENT AGAR Medium for non-fastidious bacteria growth	SUN90450	10 × 90mm	5 months
NUTRIENT BROTH Medium for the general cultivation of less fastidious microorganisms	SUN15360	Tubes: 10× 5mL	8 months
PLATE COUNT AGAR Medium for total bacterial count detection in foodstuffs and water	SUN90430	10 × 90mm	4 months
POTATO DEXTROSE AGAR Medium for the isolation and enumeration of yeasts and moulds from dairy and other food products	SUN90510	10 × 90mm	3 months
PSEUDOMONAS CN AGAR Medium for the enhancement of fluorescein production by <i>Pseudomonas spp.</i>	SUN60940	10 × 55mm	4 months
SABOURAUD DEXTROSE AGAR Selective medium for fungi and yeasts	SUN90201	10 × 90mm	4 months

SLANETZ & BARTLEY AGAR Selective medium for enterococcus isolation and enumeration	SUN60080	10 × 55mm	4 months
TRIPLE SUGAR IRON AGAR Differential medium for the identification of enterobacteria	SUN15700	Tubes: 10 × 10mL	8 months
TRYPTIC SOY AGAR Medium for the growth of a wide variety of microorganisms	SUN90260	10 × 90mm	4 months
TRYPTONE SOY BROTH Broth for the growth of a wide variety of microorganisms	SUN15340 SUN70340 SUN75340	Tubes: 10 × 5mL Bottle: 100mL Bottle: 225mL	8 months 8 months 8 months
TRYPTONE BILE X-GLUCURONIC AGAR Selective chromogenic medium for <i>E. coli</i> detection in food	SUN90490	10 × 90mm	4 months
VIBRIO (T.C.B.S) AGAR Selective medium for <i>Vibrio spp.</i> isolation	SUN90360	10 × 90mm	4 months
VIOLET RED BILE D(+)-GLUCOSE AGAR Medium recommended for the detection and enumeration of Enterobacteriaceae from food and animal feed	SUN90660	10 × 90mm	4 months
VIOLET RED BILE GLUCOSE AGAR Medium recommended for the detection and enumeration of Enterobacteriaceae from food and animal feed	SUN90320	10 × 90mm	4 months
XYLOSE LYSINE DEOXYCHOLATE AGAR (X.L.D.) Differential and selective medium for pathogenic enterobacteria isolation	SUN90280	10 × 90mm	4 months
YEAST EXTRACT AGAR Medium for microorganisms' enumeration in water	SUN15220	Tubes: 10 × 15mL	8 months
YGC AGAR Yeast Extract Glucose Chloramphenicol Agar is a ready-to-use nutrient medium specially enriched for the cultivation of yeasts and fungi in milk and dairy products	SUN90670	10 x 90mm	4 months
ORANGE SERUM AGAR Orange Serum Agar is recommended by APHA for the cultivation of Lactobacilli and other acidic organisms for citrus concentrate juices	SUN90150	10 x 90mm	3 months
YSG AGAR (Yeast Starch Glucose Agar) Yeast Starch Glucose is recommended for the enumeration of Alicyclobacillus in fruit juices	SUN90680	10 x 90mm	4 months
BRAIN HEART INFUSION AGAR Brain Heart Infusion Agar is a general-purpose nutritional medium for the cultivation of demanding, aerobic and anaerobic microorganisms including yeasts and fungi	SUN90030	10 x 90mm	4 months
GLUCOSE OF MEDIUM Glucose of Medium is a nutrient used to identify Gram (-) bacteria in products intended for human consumption, in animal feed and in environmental samples	SUN15610	Tubes: 10 x 10mL	8 months
KING'S B MEDIUM King's B Medium is a solid nutrient suitable for the detection and enumeration of Pseudomonas in water samples and is defined as the basis of the method specified in ISO 16266.	SUN15750	Tubes: 10 x 5mL	8 months
TRYPTOPHAN BROTH Tryptophan Broth is a nutrient suitable for the differentiation of indole-producing microorganisms	SUN15920	Tubes: 10 x 8mL	8 months

ΕΔΡΑ ΕΤΑΙΡΕΙΑΣ & ΣΤΟΙΧΕΙΑ ΕΠΙΚΟΙΝΩΝΙΑΣ

Sundial Co. Ltd
Χρυσάνθου Επισκόπου Πάφου 2,
3022 Λεμεσός

P.O.Box 56281
3305 Λεμεσός
Τηλέφωνο: 25 828828
Email: sundial@cytanet.com.cy
Αρ. ΦΠΑ: 100240831

ΣΗΜΑΝΤΙΚΗ ΣΗΜΕΙΩΣΗ: ΟΙ ΠΑΡΑΔΟΣΕΙΣ ΤΩΝ ΠΡΟΙΟΝΤΩΝ ΠΡΑΓΜΑΤΟΠΟΙΟΥΝΤΑΙ ΜΟΝΟ ΕΝΤΟΣ ΚΥΠΡΟΥ